

A guide to the reports

Daily tracking report

Organizes the information a member logged, including medications, symptoms, ON/OFF observations, hydration, and other recorded areas. Check the reporting period and number of recorded days before reading a summary.

Sense Suite report

Shows module measurements, baseline averages, latest results, session counts, baseline status, and practice information where recorded. VoiceSense, MotorSense, MoveSense, and MindSense cover different tasks.

Personal baseline

A baseline summarizes that member's initial recorded sessions. It is not a population norm or diagnostic threshold. Where baseline collection is incomplete, treat comparisons as unavailable.

Change and context

A numerical change describes a difference in recorded values. It does not, by itself, establish a clinical improvement or decline. Consider recording conditions, task completion, session count, and the member's observations.

A useful appointment discussion

Ask what the member noticed, which records relate to their questions, and what they should record before the next review. Patients can email or print their reports; these resources do not provide clinic access to patient accounts.

Resource status

Sense Suite practice programs are currently being vetted. These starter kits introduce tracking and report sharing; they do not supply exercise routines or claim clinical approval.